

Department	Standard Components	Room Tag:	3690
Room Name	Radiotherapy Treatment Planning	Room Code:	RAD-TRP-SJ
Net Room Area Briefed	56.0 M2		
No of Rooms	1		
Subtotal of Briefed Area	56.00 M2		
Hours of Operation	Available 24 Hours		
Occupancy	6 workstations		
Also covers these Rooms			
Description/ Special Requirements	The Treatment Planning/ Dosimetry room houses workstations and computer equipment used by the dosimetrists to determine treatment plans for each patient.		
	Acoustics: Sound attenuation - medium; Acoustic privacy is required		
Amendment	13.03.2014 - First Issue		

Room Fabric and outline specifications

ID	Fabric	Material	Finish	Specification	Selection	Remarks
5060	Ceiling	Acoustic	Prefinished	Drop-in tiles, 600x1200		
6025	Cornice	Aluminium	Powdercoat	shadow line		
8000	Door protection	Vinyl	Prefinished	Kickplate, 300mm high		
9002	Doors	Solid Core/ Glass	Paint	Single leaf, half glazed		920mm clear opening, glazing optional
1149	Floor finish	Carpet	90/10 Wool/ Nylon	Broadloom, 48 oz.		vinyl also acceptable
949496	Skirting	Timber	Paint	150mmH x 20mm thick, nominal		or floor vinyl covered skirting if vinyl used
3000	Wall finish	Paint	Acrylic, washable			

Fittings and Furniture (FF)

ID	Description	Category	Group	Qty	Selection / Remarks
1050	Bench: laminate	Furniture/ Fitting	1	2	
1120	Bench: workstation	Furniture/ Fitting	1	7	1 for server
2100	Chair: ergonomic, office	Furniture/ Fitting	3	6	
3350	Clock: battery operated	Furniture/ Fitting	3	1	
4410	Cupboard: under bench	Furniture/ Fitting	1	1	
4460	Cupboard: wall mounted	Furniture/ Fitting	1	1	over bench
17255	Pedestal filing unit: mobile, 660mmH	Furniture/ Fitting	3	6	under each workstation
17360	Pinboard: fabric covered	Furniture/ Fitting	1	6	above workstation benches
20000	Shelving unit: adjustable	Furniture/ Fitting	1	6	above workstation benches
24100	Waste bin: paper, 20 Litre	Furniture/ Fitting	3	6	under workstations
24400	Whiteboard	Furniture/ Fitting	2	1	900x900, optional

Fixtures, Equipment and associated Services (FE)

ID	Description	Group	Qty	Ele	Data	Cd	WH	TW	/m	WT	Tap	Dns	Gas	Selection / Remarks
														
Provided by TAHPI for : Sharjah Domain					Issue No: 1 Issue Date: 1 Apr 2014					Print Date: 20 May 2014 Page 1 of 2				

Department	Standard Components	Room Tag:	3690
Room Name	Radiotherapy Treatment Planning	Room Code:	RAD-TRP-SJ

11400	Computer: PC & monitor	3	6	yes	yes	___	___	___	___	___	___	___
11900	Computer: server	3	1	yes	yes	___	___	___	___	___	___	___
2039	Film digitiser	3	1	yes	yes	___	___	___	___	___	___	___
41100	Printer: networked	3	1	yes	yes	___	___	___	___	___	___	___
54510	Telephone: handset, standard	3	6	___	yes	___	___	___	___	___	___	___

Services

ID	Description	Service Category	Qty	Selection / Remarks
1006	Voice / Data outlet: double	Communications	6	1 outlet each workstation for telephone
1006	Voice / Data outlet: double	Communications	3	1 for printer; 1 to server; 1 on bench
5000	Airconditioning	HVAC		
6000	General fluorescent	Lighting		
6055	Special: task light, built in	Lighting		under O/H joinery
9001	GPO: Double	Power	16	2 per workstation + 1- printer
9060	GPO: Special: for equipment	Power	2	for XRV