

<b>Department</b>	<b>Standard Components</b>	<b>Room Tag:</b>	<b>2100</b>
<b>Room Name</b>	<b>CT Scanning - Control Room</b>	<b>Room Code:</b>	<b>CTCR-SJ</b>
<b>Net Room Area Briefed</b>	6.0 M2		
<b>No of Rooms</b>	1		
<b>Subtotal of Briefed Area</b>	6.00 M2		
<b>Hours of Operation</b>	Available 24 Hours		
<b>Occupancy</b>	Up to 2 staff		
<b>Also covers these Rooms</b>			
<b>Description/ Special Requirements</b>	Room adjacent to CT Scanning Room/s with direct visibility to the patient for remote control of equipment, and review and reporting of procedure images. The Control room may serve two procedure rooms.  Acoustics: Acoustic privacy is required  Additional Design Considerations: - Radiation shielding as advised by Radiation Consultant - Cable ducts from CT Scanning and Computer room may be required		
<b>Amendment</b>	13.03.2014 First issue		

**Room Fabric and outline specifications**

ID	Fabric	Material	Finish	Specification	Selection	Remarks
5060	Ceiling	Acoustic	Prefinished	Drop-in tiles, 600x1200		
6025	Cornice	Aluminium	Powdercoat	shadow line		
8000	Door protection	Vinyl	Prefinished	Kickplate, 300mm high		
9004	Doors	Solid Core/ Glass	Paint	Single leaf, observation panel		staff entry/ exit; 920mm clear opening
1000	Floor finish	Vinyl	Standard slip resistant	Seamless, coved		
7025	Observation	Lead glass	Clear	Internal window		to procedure room
2000	Skirting	Vinyl	Prefinished	Floor vinyl coved, 150 mmH		
3000	Wall finish	Paint	Acrylic, washable			

**Fittings and Furniture (FF)**

ID	Description	Category	Group	Qty	Selection / Remarks
2100	Chair: ergonomic, office	Furniture/ Fitting	3	2	
3400	Coathook	Metalwork	1	2	
22615	Table: computer	Furniture/ Fitting	3	1	small
24400	Whiteboard	Furniture/ Fitting	1	1	1200 x 900 optional

**Fixtures, Equipment and associated Services (FE)**

ID	Description	Group	Qty	Ele	Data	Cd	W	H	T	W	m	W	Tap	Dns	Gas	Selection / Remarks
11400	Computer: PC & monitor	3	1	yes	yes	___	___	___	___	___	___	___	___	___	___	
62150	IMG: console, control	1	1	yes	yes	___	___	___	___	___	___	___	___	___	___	CT Scanning equipment

Department | Standard Components | Room Tag: 2100  
 Room Name | CT Scanning - Control Room | Room Code: CTCR-SJ

ID	Description	Group	Qty	Ele	Data	Cd	WH	W	m	WT	Tap	Dns	Gas	Selection / Remarks
54510	Telephone: handset, standard	3	1		yes									

Services

ID	Description	Service Category	Qty	Selection / Remarks
1040	Intercom	Communications	1	with handsfree microphone to CT Scan room
1006	Voice / Data outlet: double	Communications	3	1 outlet for telephone
5000	Airconditioning	HVAC		
6000	General fluorescent	Lighting		
9001	GPO: Double	Power	5	
9060	GPO: Special: for equipment	Power	1	CT Scanning equipment
12030	Lead lining	Shielding		to wall shared with CT Procedure Room